# ENDO STORE

Endo & More




# ABOUT US!

Our Company, where dental care meets excellence. At Endo & more, we take pride in providing top-notch dental supplies and marketing solutions.

**Our Mission** is to support dental professionals in delivering exceptional patient care by offering a wide range of premium products and effective marketing strategies.

Join US on this journey as we strive to elevate the standard of dental practices worldwide







minfo@endostore.ae







ENDO STORE AGENTS









**Endo With Eighteeth** 



E-FLOW

DENTAL ANESTHESIA

DERLIVERY SYSTEM





## HELIOS 500





## Why HELIOS 500?

#### Fit for your budget

The price of HELIOS 500 is affordable for almost everyone, and at the same time, it provides high performance.

#### Improve your digital workflow

One of the most accurate and fastest scanners on the market.

#### Free your hands

Automatic and adaptive calibration, manual calibration is not necessary any more.

#### Make your experience comfortable

Smaller and lighter while achieving a 16mm x 14mm high-definition field of view.

More friendly to people with sensitive teeth.

(The air duct inside the body has been upgraded.)



## HELIOS





### Accuracy

As one of the most accurate scanners on the market, HELIOS 500 has reached the state-of-the-art accuracy level- the average deviation of a in-vivo full arch scan is around 20.0 µm.




Speed
Scan a full arch in less than 30 seconds with HELIOS 500.



## Weight

Only 185 grams, one of the most lightweight scanners on the market, well-balanced, ergonomically designed, and easy to carry.





## Al technology

Al technology filters the soft and moving tissues, offering users a smooth scanning experience and generating a natural-looking treatment outcomes.

# Patients Are Willing To Pay For Less Pain!




DENTAL ANESTHESIA DERLIVERY SYSTEM



## High-precision Liquid Control

While achieving precise drug delivery accuracy, you can also see the real-time remaining dose

#### Real-time Dynamic Pressure



Reflecting the injection pressure based on the value and dynamic bar graph, of which dentists can adjust the operation to bring a painless and safe treatment





#### **Thoughtful Training Mode**

Provide additional voice reminders not available in the standard mode, help users quickly learn the proper operation



#### **Smart Voice Reminder**

Timely voice reminders under low anesthesia dose and high system pressure (in PDL mode)

## 4 injection modes,

#### multiple options for effective pain relief













#### PDL Mode (periodontal ligament)

Automatic injection can be selected in this mode

#### **Dosage Mode**

Free choice of injection dose, making the injection more stable and precise









## PREMIXED INJECTABLE

## ENDOCEM MTA PREMIXED REGULAR



[Product composition]
2g Syringe 1EA / 22gauge Tip 10EA / 20gauge Tip 10EA

The first premixed MTA based on pure C3S with excellent physical properties and biocompatibility comparable to powder type MTA

Endocem MTA Premixed Regular

PERFORATION

**APEXOGENESIS** 

**PULPOTOMY** 

- Hemostatic action by phyllosilicate ingredient.
- Easy to use without bleeding control.
- High compressive strength and stable dimensional change.
- Aesthetic by Zirconium Oxide.
- Easy application to narrow and deep areas.
- Outstanding biocompatibility



## ENDOSEAL MTA



ENDOSEAL MTA is a paste-type root canal sealer based on pozzolan cement that has excellent physical and biological properties of MTA. It is premixed and pre-loaded in a syringe that allows direct application of the sealer into the root canal without requiring powder/liquid mixing, unlike other products. The product has outstanding flowability and maneuverability, which makes it possible to completely fill the root canal system including accessory and lateral canals. It is eugenol-free and will not impede adhesion inside the root canal.

**ENDOSEAL** MTA is developed for root canal filling and root perforation repair, effectively preventing periapical tissue irritation as well as secondary infection.

## STATE-OF-THE-ART MTAROOT CANAL SEALER THAT OFFERS

- Fast setting time, 10-12 minutes
- Strong antibacterial effects
- Superior dimensional stability
- Easy retrieval with NiTi files
- Hermetic apical sealing
- High radiopacity
- Promotion of hard tissue formation
- No postoperative pain after overfilling





## ONES MED






**IZY FIL** 



**IZY FLOW** 



**IZY ETCH** 



**IZY CORE** 



**IZY CEM** 



**IZY PREP** 



**IZY TEMP UNIVERSAL** 



**IZY TEMP CROWN** 



## Izy Fil

**Light Cured Composite Resin** 

Packing 4g x 1ea
Shade A1, A2, A3

**RESTORATIVE MATERIALS** 

## **Izy Flow**

Light Cured Flowable Resin

Packing 1.2ml(2g) x 2ea

Shade A2





## Izy Etch

**Etching Agent** 

Packing | 5ml x 2ea or 4ea, Tips



#### **RESIN CEMENT MATERIALS**

## **Izy Core**

dual cured core build up resin

Packing | 5ml x 2ea, Tips





## Izy CEM Universal

Dual Cured Permanent Resin Cement

Packing | 5ml(8g) x 2ea

## Izy Prep

**EDTA Cream Cleaner** 

Packing | 5ml(6g) x 2ea, Tips

**ENDODONTIC MATERIALS** 



#### **TEMPORARY MATERIALS**

## Izy Temp Universal

Light Cured Temporary Filling Material

Packing 3g x 2ea

Shade Tooth color(A2)





Self Cured Temporary Crown & Bridge

Packing 50ml x 1ea
Shade A1, A2, A3

**Izy Temp Crown** 



Quality
Product Development



Suzhou RAMO Medical Instrument Co., Ltd

**RCS RAINBOW ONE FILE** 

RCS MAGIC PATH ROTARY FILE

RCS GOLD LT ROTARY FILE

RCS RETREATMENT ROTARY FILE

RCS CHILDREN FILE

**RCS SHAPE D** 

RCS K FILE



## RCS RAINBOW

ONE FILE One file system, simple and efficient instrument



**FEATURES** 

- Flexible
- Upgraded slim S-shaped cross section.
- Special heat-treated Tech and surface treatment, reduce crack and improve resistance to cyclic fatigue.
- Greater torsional resistance and flexural strength.
- Improved cutting efficiency by more than 150%.

# RCS MAGIC PATH ROTARY FILE Simple | Safe | Efficient



- Rotary path file.
- **FEATURES**
- Suitable for curved canals.
- More flexible than K-files.
- Faster instrument time than manual K-files.
- 7 sizes available.

SPEED&TORQUE 350RPM/1.5N



# RCS GOLD LT ROTARY FILE

Simple | Safe | Efficient

FEATURES

- Pro-Taper system.
- OSX, S1,S2 are flexible; F1, F2, F3 are extremely flexible,control memory.
- Make sure shaping file have good cutting efficiency and finishing file can save more dentin.



# RCS RETREATMENT ROTARY FILE Simple | Safe | Efficient



#### **FEATURES**

High shaping to remove gutta-percha and paste fillers fast.
5 sizes meet all retreatment request.

SPEED&TORQUE 350RPM/2.0N



## RCS CHILDREN FILE



#### **FEATURES**

- Super flexible.
- Short length design is better for children use.
- Improved cutting efficiency and resistance of cyclic fatigue.

## RCS SHAPE D

Simple | Safe | Efficient



- Stainless steel, rigid
- Shaped D cross section
- Use for enlarging root canal and negotiating calcified and narrowed root canal

### **FEATURES**





